



Are Sparging Stones Recommended for Ozonating Water - FAQ

Date: 6-APRIL-2012 Last Updated: 8-NOVEMBER-2025

Yes, you can use any of our stainless steel sparging stones with ozone to "ozonize" water. It is recommended to use 10um pore size stones.

[Click here for ordering information and pictures of our sparging stones](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com