

What is the pressure limit of the stainless steel filter sparging stones - FAQ

All the sparging stone products are rated to a maximum pressure of 1,000 bar used under standard conditions.

However, it is important to note that the tubing used and the connection of the tubing to the device may not perform at this pressure rating. The limiting pressure factor of a filter sparging stone is not for the filters / spargers, but for the connections.

Click [HERE](#) for filter sparging stone ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com