



What is the pH Range of 316 Stainless Steel Tubing Fittings and Sparging Stones - FAQ

Date: 12-MAY-2012 Last Updated: 3-NOVEMBER-2025

The useful pH range of 316 stainless steel tubing, fittings and sparging stones is 1 - 14.

Click [HERE](#) for more information about stainless steel tubing

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com