



## Are Stainless Steel Tubing or Fittings Autoclavable - FAQ

*Date: 12-MAY-2012 Last Updated: 3-NOVEMBER-2025*

Yes, many scientists have successfully used 316 stainless steel after autoclaving tubing and or fittings with gamma irradiation, heat or ethylene oxide sterilization methods.

Click [HERE](#) for stainless steel tubing ordering information

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)