

If I am using TFA in HPLC will I have a problem with my Delrin type fittings - FAQ

Delrin is not compatible with trifluoroacetic acid (TFA).

If it comes in contact with your fittings, your fitting could become compromised and your connection might slip and become less than perfect or acceptable.

Click [HERE](#) to view our HPLC fittings

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com