

Choosing the correct HPLC column coupler - How To

To find the correct size of the column coupler, simply match your system tubing's internal dimension (ID) to the column coupler's ID.



Click [HERE](#) for HPLC column coupler ordering information

The different coupler colors correspond with those of the HPLC tubing products .

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com