

What is the Definition of a Fitting for HPLC - FAQ

Date: 7-APRIL-2012 Last Updated: 8-NOVEMBER-2025

The dictionary defines a fitting as “a small, often standardized part used to join other parts together properly.”

This is very similar to what plumbers refer as fittings. For the use of chromatographers, fittings consist of **nuts with internal and external threads** for both high pressure and low pressure connections. Also, unions, crosses and tees and adapters are considered fittings for both high and low pressure connections.

Click [HERE](#) to view our HPLC fittings

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com