

What Is the Useful Temperature Range of PEEK Laboratory Products - FAQ

Date: 12-MAY-2012 Last Updated: 28-NOVEMBER-2025

The useful temperature range of PEEK related products for chromatography is up to 100°C (tubing) and up to 125°C (fittings).

PEEK is sometimes rated to 250°C but... It will not be as stable (lower burst pressure for tubing) as it is normally and some solvents may interact stronger than under RT with PEEK.

NOTE: PEEK has to be used with caution with higher temperatures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com