



What is the Useful Temperature Range of 316 Stainless Steel Laboratory Products - FAQ

Date: 12-MAY-2012 Last Updated: 8-NOVEMBER-2025

The useful temperature range of stainless steel (SS) related products such as column hardware, tubing and HPLC fittings is up to 240°C.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com