



What is the pressure rating of the Hichrom brand column coupler - FAQ

The Hichrom brand column coupler catalog number 47450-06 is rated for use up to 6,000 psi.

Click [HERE](#) for column coupler ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com