



## Does the Direction of a Guard Column in the Holder Matter - FAQ

**No, the guard column cartridge can be placed in the holder in either direction.**

The label on the column does not indicate a particular direction to place the column in the holder. The guard columns have no flow direction in the holder.

[Click HERE for Ordering Information Including Guard Column Holders](#)



Printed from the Chrom Resource Center  
Copyright 2025, All Rights Apply  
**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900  
Fax: (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)