



Stability of Silica-C Columns Compared to Other Silica Columns - Tech Information

We have done extensive testing in our labs and in real world use around the world and have proven that the Cogent Silica-C™ columns are stable even in the presence of water due to the silica hydride's slightly hydrophobic surface.

However, these columns are designed and specified for Normal Phase and Aqueous Normal Phase ANP HPLC and not for Reversed Phase. But the only real limitation is your creativity.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com