



Are the Cogent TYPE-C Silica Stationary Phases Available for Sale as Bulk Materials - FAQ

Date: 16-AUGUST-2012 Last Updated: 8-NOVEMBER-2025

No, the Cogent 4um and 2.0 TYPE-C Silica™ based stationary phases are not available for sale as unpacked bulk material.

TYPE-C Silica™ is a very unique material and to properly pack columns or capillaries requires proprietary technology. For this reason, it is only sold in a pre-packed column format. However, 10um particles can be made available upon request.



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com