



Do You Offer a 0.5um Titanium or PEEK Frit for the PEEK Column Hardware - FAQ

Date: 2-MAY-2013 Last Updated: 28-NOVEMBER-2025

In addition to the standard 5um Stainless Steel or PEEK frits, we have 0.5um Titanium Frits available, but not in PEEK .

[PEEK Column Hardware Product Page](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com