

## Does the Hichrom brand column coupler work with the Partisil and Partisphere HPLC columns - DISCO

*Date: 19-APRIL-2017 Last Updated: 14-NOVEMBER-2025*

**Yes!**

The column coupler with the catalog number [47450-06](#) is a suitable connector for the Partisil® & Partisphere® brands of HPLC columns or to other brands of columns or to a guard column with compatible threads.



Column Coupler for use with guard or analytical columns.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)