

## Can Chemical Solutions Qualification Kits Be Used with Any HPLC Instrument - FAQ

*Date: 27-SEPTEMBER-2012   Latest Update: 23-JUNE-2025*

✔ **Yes!** MICROSOLV's Chemical Solutions brand HSQ™ and PQ™ Kits are compatible with ANY HPLC System.

Looking for a flexible, high-performance qualification solution for your HPLC system? The MICROSOLV HSQ™ and PQ™ Performance Qualification Kits are designed to work seamlessly with any brand of HPLC instrument and any chromatography data system — including Waters Empower™.

### **Empower™ Users: What You Should Know**

While Waters Corporation offers a PQ kit that integrates directly with Empower™ for automated execution, the MICROSOLV kits offer a broader, more comprehensive qualification. Our kits are not controlled by Empower™, so they won't run automatically within that software. However, you can still:

- Use Empower™ to acquire and integrate chromatographic data
- Easily transfer results into the included Excel-based software provided with our kits for analysis and documentation

### **Why Choose Chemical Solutions Qualification Kits?**

- ✔ Broader qualification scope than Waters' kit
- ✔ NIST-traceable materials for confidence and compliance
- ✔ Compatible with all HPLC systems and data platforms
- ✔ Detailed, Excel-based reporting tools included
- ✔ Trusted by labs worldwide for over 25 years

Whether you're qualifying a single instrument or managing a fleet across multiple labs, MICROSOLV's Chemical Solutions Qualification kits provide the flexibility, accuracy, and depth you need — without locking you into a single software ecosystem.

---

### ***Range of Benefits of the Chemical Solutions Qualification Kits.***

#### ✔ **1. Universal Compatibility**

- Works with any brand of HPLC or UV-Vis instrument
- Compatible with all chromatography data systems, including Waters Empower™, Agilent OpenLab, Chromeleon, and more

## 2. NIST-Traceable Standards

- All solutions are NIST traceable, ensuring accuracy, reliability, and regulatory confidence
- Ideal for labs that require audit-ready documentation

## 3. Comprehensive Qualification Scope

- More extensive than many OEM kits, covering:
  - Detector linearity
  - Wavelength accuracy
  - Flow rate precision
  - Injection volume accuracy
  - System suitability and more
  - Supports compliance with USP <1058> and USP <857>

## 4. High-Quality, Ready-to-Use Solutions

- Pre-measured, pre-labeled, and lot-specific
- Long shelf life and consistent performance
- Reduces preparation time and human error

## 5. Expertly Designed by Industry Leaders

- Developed with input from Dr. Richard Hartwick, formerly of Chemical Solutions, Inc.
- Backed by MICRO SOLV's technical support team with decades of experience

## 6. Excel-Based Software Included

- Easy-to-use spreadsheets for data entry, analysis, and documentation
- No proprietary software required
- Compatible with data exported from any CDS

## 7. Cost-Effective & Flexible

- No need for expensive software licenses or locked-in systems
- Use only the components you need — ideal for labs with multiple instrument types

 Click [HERE](#) for PQ & HSQ Kit ordering information and product images

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)