

Qualification Kits Use with Most HPLC Instruments - Tech Information

Date: 27-SEPTEMBER-2012 Latest Update: 23-JUNE-2025

✔ MICRO SOLV's Chemical Solutions brand HSQ™ and PQ™ Kits are compatible with ANY HPLC System.

Looking for a flexible, high-performance qualification solution for your HPLC system? The MICRO SOLV HSQ™ and PQ™ Performance Qualification Kits are designed to work seamlessly with any brand of HPLC instrument and any chromatography data system — including Waters Empower™.

Empower™ Users: What You Should Know

While Waters Corporation offers a PQ kit that integrates directly with Empower™ for automated execution, the MICRO SOLV kits offer a broader, more comprehensive qualification. Our kits are not controlled by Empower™, so they won't run automatically within that software. However, you can still:

- Use Empower™ to acquire and integrate chromatographic data
- Easily transfer results into the included Excel-based software provided with our kits for analysis and documentation

Why Choose Chemical Solutions Qualification Kits?

- ✔ Broader qualification scope than Waters' kit
- ✔ NIST-traceable materials for confidence and compliance
- ✔ Compatible with all HPLC systems and data platforms
- ✔ Detailed, Excel-based reporting tools included
- ✔ Trusted by labs worldwide for over 25 years

Whether you're qualifying a single instrument or managing a fleet across multiple labs, MICRO SOLV's Chemical Solutions Qualification kits provide the flexibility, accuracy, and depth you need — without locking you into a single software ecosystem.

Range of Benefits of the Chemical Solutions Qualification Kits.

✔ **1. Universal Compatibility**

- Works with any brand of HPLC or UV-Vis instrument
- Compatible with all chromatography data systems, including Waters Empower™, Agilent OpenLab, Chromeleon, and more

2. NIST-Traceable Standards

- All solutions are NIST traceable, ensuring accuracy, reliability, and regulatory confidence
- Ideal for labs that require audit-ready documentation

3. Comprehensive Qualification Scope

- More extensive than many OEM kits, covering:
 - Detector linearity
 - Wavelength accuracy
 - Flow rate precision
 - Injection volume accuracy
 - System suitability and more
 - Supports compliance with USP <1058> and USP <857>

4. High-Quality, Ready-to-Use Solutions

- Pre-measured, pre-labeled, and lot-specific
- Long shelf life and consistent performance
- Reduces preparation time and human error

5. Expertly Designed by Industry Leaders

- Developed with input from Dr. Richard Hartwick, formerly of Chemical Solutions, Inc.
- Backed by MICRO SOLV's technical support team with decades of experience

6. Excel-Based Software Included

- Easy-to-use spreadsheets for data entry, analysis, and documentation
- No proprietary software required
- Compatible with data exported from any CDS

7. Cost-Effective & Flexible

- No need for expensive software licenses or locked-in systems
- Use only the components you need — ideal for labs with multiple instrument types

 Click [HERE](#) for PQ & HSQ Kit ordering information and product images