

Shorter Column Supplied in the HSQ Kit Benefits - Tech Information

Date: 9-APRIL-2014 Last Updated: 1-JULY-2025

The HSQ Kit™ comes with a 4.6 x 75mm column that is used in the qualification procedure.

The older version of the the HSQ Kit™ (called the PQ Kit™) comes with a 4.6 x 150mm column instead.

This difference will not affect the qualification data results and the columns can be used interchangeably. The only difference for the end-user is that the HPLC runs take much less time due to the shorter column length when using the HSQ Kit™. Hence, this anadvantage over the older PQ Kit.

Click [HERE](#) for HSQ & PQ kits ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com