



Storing the Chemical Solutions Qualification Kit Solutions When Not in Use - How To

Date: 7-APRIL-2012 Last Updated: 2-JULY-2025

The test solutions should be stored at room temperature **20-25°C** in a safe, dry place. Never refrigerate these solutions.

Click [HERE](#) for Chemical Solutions PQ & HSQ Qualification Kit ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com