



Storing the HPLC column supplied in the Chemical Solutions Qualification kits - How To

Date: 7-APRIL-2012 Last Updated: 2-JULY-2025

One of the best features of the Chemical Solutions Qualification Kits™ is that these DIY kits were designed for ease and speed. The certified, special HPLC column is shipped and stored in the same mobile phase that is used during the testing.

The column is always ready to go with minimal equilibration time. The mobile phase / storage solvents were chosen for long column life. The solutions should be stored at **room temperature** in a dry place and never refrigerated.

Click [HERE](#) for Chemical Solutions PQ Kit™ and HSQ Kit™ Ordering Information and Pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com