

How Can I Obtain Qualification Solutions with a Longer Expiration Date Than Normal - FAQ

Date: 13-SEPTEMBER-2013 Last Updated: 30-JUNE-2025

All of our Chemical Solutions brand of Qualification Solutions are manufactured with a standard two-year expiration date, supported by validated stability data. Both the date of manufacture and the expiration date are clearly printed on the Certificate of Analysis (COA) included with each product.

You can use these solutions with full confidence up to the expiration date, and even beyond, if your internal SOPs allow for it.



Need a Longer Expiration Date?

If you require qualification solutions with a longer shelf life than what is typically shipped, we can provide a custom-made batch. Please note:

- Lead time: Approximately 4- 5 weeks ARO
- Minimum order quantity (MOQ) may apply



Shipping Policy

As a standard practice, MICROSOLV ships products with at least six months of shelf life remaining, ensuring most customers have ample time to use the solutions before expiration.



Click [HERE](#) for HSQ & PQ Kits ordering information and product images.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com