

How many tests will I get out of one of the Chemical Solutions Qualification Kit - FAQ

Date: 11-APRIL-2012 Last Updated: 28-JUNE-2025

Each of the Chemical Solutions HPLC Qualification Kits™ are designed to support up to 30 qualification tests, depending on how carefully the solutions are handled and transferred. These tests can be performed on:

- A single instrument multiple times, or
- Up to 30 different instruments, using one test per system.

💡 Helpful Tip:

To maximize the number of tests:

- Use precise pipetting techniques
- Minimize waste during solution transfer
- Store reagents according to the kit instructions

 Click [HERE](#) for PQ Kit™ and HSQ Kit™ ordering information and product images

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com