

Catalog Number

45208-F-ETFE

Ferrules, for use with low pressure, tubing connectors. Use with flange free, 1/8 inch, low pressure tubing. ETFE. ARE-Applied Research brand. 10/PK.



Additional Info

Use with tubing connector nuts for low pressure tubing connections in HPLC or other fluidic connections.

ETFE (Ethylene tetrafluoroethylene) is a fluorine-based plastic. It was designed to have high corrosion resistance and strength over a wide temperature range. ETFE has a relatively high melting temperature, excellent chemical resistance properties.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$26.00

Printed On 2024/05/20