

## Catalog Number

07210-M

**Magnifier, 15x, Working Distance 0.5", MicroSolv Brand. 1 EA.**



### Additional Info

This 15x Magnifier is often used to Inspect CE or CEC Capillaries for air bubbles or cracks in the Fused Silica that can cause poor Data or loss of Current.

The working distance is 0.5" and the unit is only 1.375" H x 0.75" W at the base with a 0.5" Field of View

The transparent Base allows light in so lamps may not be necessary but when using a light box, you can see directly into the Capillaries.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$196.86**

Printed On 2024/04/30