

Catalog Number

47001-002

Sample Loop, Cheminert/Rheodyne, 2 ul, For HPLC Valves, PEEK.
ARE-Applied Research brand. 1 EA.



Additional Info

These HPLC injector sample loops are made with PEEK Tubing for unsurpassed injection reproducibility and accuracy.

Tubing is 0.06mm ID (+/- 0.025mm) and is supplied with 2 One-Piece, Small Hex Head, PEEK, High Pressure Tubing Connectors for fitting on the injectors. A label with the loop volume is attached.

Use these Loops with Rheodyne® or Cheminert® Brand Injector Valves.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$46.23

Printed On 2024/02/25