

Catalog Number

47004-250

Sample Loop, Valco, 250 μ l, For GC Injectors, Stainless Steel. ARE-
Applied Research brand. 1 EA.



Additional Info

These injector sample loops are made with welded and annealed, 316 stainless steel tubing for unsurpassed injection reproducibility and accuracy.

Tubing is 0.75mm ID (+/- 0.025mm) and is supplied with 2 two-piece, small hex head, SS, high pressure tubing connectors with ferrules for fitting on the injectors. Tubing end cuts are made using EDC (electronic discharge) techniques to achieve square and burr free ends for Zero-Dead Volume Connections. A label with the loop volume is attached.

Use these loops with VICI® or Valco® brand GC injector valves.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$59.70

Printed On 2024/05/02