

Catalog Number

49302-CX-09

Union, Connector Kit, Cross, high pressure, 0.8 mm / 0.031 inch through hole with 1/4-28 screw threads. PEEK, round with an internal "X" flow pattern. Use with 1/16 inch tubing. ARE-Applied Research brand.



Additional Info

This "Union" is used to connect PEEK or stainless steel tubing the same connectors on all ends in the high pressures zones of HPLC systems. A "Cross (X)" internal flow pattern is often used to split flow into two or three directions or to combine two or three streams into one.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/06