

Catalog Number

49302-CX-09

Union, Connector Kit, Cross, high pressure, 0.8 mm / 0.031 inch through hole with 1/4-28 screw threads. PEEK, round with an internal "X" flow pattern. Use with 1/16 inch tubing. ARE-Applied Research brand.



Additional Info

This "Union" is used to connect PEEK or stainless steel tubing the same connectors on all ends in the high pressures zones of HPLC systems. A "Cross (X)" internal flow pattern is often used to split flow into two or three directions or to combine two or three streams into one.

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com