

Catalog Number

49302-CX-11-005

Union, Connector Kit, Cross, high pressure, 0.30 mm / 0.012 inch through hole with 10-32 screw threads. PEEK, round with an internal "X" flow pattern. Use with 1/16 inch tubing. ARE-Applied Research brand.



Additional Info

This low dead volume, "Union" is used to connect PEEK or stainless steel tubing with the same connectors on both ends in the high pressure zones of HPLC systems. A "Cross (X)" internal flow pattern is often used to split flow into two or three directions or to combine two or three streams into one.

Supplied with three finger tighten connectors (nuts) [47301-05](#).

Note: Click [HERE](#) for key differences between high and low pressure unions and adapters.

Price

\$811.61

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$811.61