

## Catalog Number

49302-T10-010

**Union, Connector, Tee**, high pressure, 0.30 mm / 0.012 inch through hole with 10-32 screw threads. PEEK, round with an internal "T". ARE-Applied Research brand. 10/PK.



### Additional Info

This "Union" is used to connect PEEK or stainless steel tubing on all ends, in the high pressures zones of HPLC systems. The "Tee Connection" internal flow pattern is often used to split flow into two directions or to combine two streams into one.

Click [HERE](#) for suggested finger tightened connectors (nuts) which are sold separately.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,136.21**

Printed On 2024/05/16