

## Catalog Number

49302-T11-010

**Union, Connector Set, Tee**, high pressure, 0.30 mm / 0.012 inch through hole with 10-32 screw threads. PEEK, round with an internal "T". Use with 1/16 inch tubing. ARE-Applied Research brand. 10/PK.



### Additional Info

This "Union" set is used to connect PEEK or stainless steel tubing on all ends, in the high pressures zones of HPLC systems. A "Tee Connection" is often used to split flow into two directions or to combine two streams into one.

Dead Volume - 0.689  $\mu$ L

Supplied with 3 finger tighten connectors (nuts) [47320-10N](#) .

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,001.53**

Printed On 2024/05/17