

Catalog Number

49681-05-005

Inline Filter, for HPLC, female connections replaceable filters and All-PEEK body. Includes 25 disposable 0.5 um SS filters. ARE-Applied Research brand. 5/PK.



Additional Info

Compatible with typical flow rate analytical columns up to 4,500 psi and 5 ml/minute. Also used as a union to connect tubing, but with a filter. *Not recommended for use with UHPLC pressure ranges.*

- Economical way to protect your columns and instruments.
- Re-usable filter holder.
- Easy to install after the pump, before the injector.
- Easy to replace the filter.

Price**\$237.73**

Additional Info

The universal connection (fits all standard HPLC threads and pilot lengths) can be installed with a wrench or hand tightened. When these filters are in use, and pressure starts to rise, simply disconnect (after pressure in the column is at atmospheric), dispose of and replace with a new filter.

A pack of 25 each, 0.188 inch, 0.5 um stainless steel, PEEK encased filters are supplied but other PEEK encased filters with the same or different porosity, that are 1/4 inch OD and 1/16th inch thick should work with this holder.

Click [HERE](#) for replacement filter ordering information. Click [HERE](#) for economical pack of 10 PEEK holders and 25 SS filters.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$237.73