

Catalog Number

05080-W

CE Water, Ultrapure for Capillary Electrophoresis. 80ml, MicroSolv Brand. 5/PK.



Additional Info

This 18 megohm, Highly Purified Deionized Water is for use in CE, CZE or any other use that requires purity and virtually free of Background Electrolytes (BGE).

Produced in Polypropylene Bottles.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$67.66

Printed On 2024/05/19