

Catalog Number

40000-05P

HPLC Column, Silica-C, 4um, 4.6mm ID x 50mm Length, 100A. Cogent
TYPE-C Brand. 1 EA.



Additional Info

The Silica-C™ HPLC Columns are made with [Cogent TYPE-C™ Silica](#), and provide unique Selectivity due to its surface chemistry (especially in Normal Phase HPLC) compared to standard silica, This Silica Hydride Phase is not modified and therefore has no end-capping or organic ligands.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

This size column provides good resolution and short run times and is a Very Popular Column used in High Pressure Through-put Screening.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$839.68

Printed On 2024/05/20