

Catalog Number

40000-P21-25-3M

Prep HPLC Column, Silica-C, 5um, 300A. 21.2mm ID x 250mm Length, Preparative. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Silica-C™ Preparative, HPLC Column is made with 5um, 300A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent yields and turn around times.

Chosen for analyzing large compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes
- Use Minimal Salt to Elute in ANP.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$7,234.25

Printed On 2024/05/20