

Catalog Number

40008-SP250-3M

Semi Prep HPLC Column, Bidentate C8 300, 5um, 300A. 10mm ID x 250mm Length. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Bidentate C8 300 "Semi-Preparative", HPLC Column is made with 5um, 300A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Very Robust Column Lifetimes.

Chosen for Purification or "Prep" of Hydrophobic Macro-Molecules such as Proteins and Poly Saccharidesm when a large volume of sample is desired to be collected.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$3,555.16

Printed On 2024/05/20