

## Catalog Number

40031-PG5

**Guard Column cartridges, UDA (wcx)**, replacement, 4.0 mm ID x 10 mm long cartridges packed with 4  $\mu$ m, 100 Å phase. Cogent TYPE-C Brand. 5/PK.



## Additional Info

These Guard Columns are designed to protect your analytical columns from particulate and contaminants that will adsorb to your column. These columns are packed in Stainless Steel tubes with SS Frits. Labelled, and shipped in individual, black, delrin cases to protect and keep them from drying out. Use with [guard column holder](#) sold separately.

*Use this guard column cartridge with a 3.0, 4.0, and 4.6 mm ID analytical columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$889.59**

Printed On 2024/05/20