

## Catalog Number

40036-05P-2

**HPLC Column, Amide, 4um, 2.1mm ID x 50mm Length, 100A. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Amide HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of Sugars and other Analytes with amine groups,

*This size column provides short run times but requires a pump that can precisely deliver, and a flow cell that made for, reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$899.43**

Printed On 2024/05/20