

Catalog Number

40036-10P

HPLC Column, Amide, 4um, 4.6mm ID x 100mm Length, 100A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Amide HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of Sugars and other Analytes with amine groups,

This size column is often used when shorter run times are needed but 50mm long columns do not provide enough resolution.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,064.53

Printed On 2024/05/20