

## Catalog Number

40036-P30-50

**Prep HPLC Column, Amide, 4um, 100A. 30mm ID x 50mm Length.**  
Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Cogent Amide preparative, HPLC column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent selectivity, faster equilibration and robust column lifetimes.

Chosen for purification or "Prep" of sugars and other analytes with amine groups when a large volume of sample is desired to be collected.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$15,470.56**

Printed On 2024/05/20