

Catalog Number

40100-05P-2

Method Development Kit 2054, TYPE-C. Includes 5 Cogent HPLC Columns, 2.1mm x 50mm. Cogent TYPE-C Brand. 1 PK.



Additional Info

This Kit, 2054 is an economical way to purchase 5 Cogent HPLC Columns that are often used with LCMS for Method Development and Scouting. Each Column has a 2.1mm Internal Diameter and is 50mm Long, packed with 4um particles, 100A pore size and are manufactured with [Cogent TYPE-C™ Silica Hydride](#).

This kit includes the columns listed below:

- **Cogent Diamond Hydride™**, Click [HERE](#) for More Information.
- **Cogent Bidentate C18™**, Click [HERE](#) for More Information.
- **Cogent Bidentate C8™**, Click [HERE](#) for More Information.
- **Cogent Silica-C™**, Click [HERE](#) for More Information.
- **Cogent UDC-Cholesterol™**, Click [HERE](#) for More Information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$4,150.79

Printed On 2024/05/20