

Catalog Number

40200-75P-2

HPLC Column, Silica-C 2.0, 2.2um, 2.1mm ID x 75mm Length, 120A, High Pressure. Cogent TYPE-C Brand. 1 EA. CUSTOM ORDER



Additional Info

This Cogent Silica-C™ “near UHPLC”, HPLC Column is made with 2.2um, 120A, Cogent [TYPE-C™ Silica](#), and this silica hydride based stationary phase.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require nonpolar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,163.03

Printed On 2024/05/20