

Catalog Number

40218-HG1

Guard Column Cartridges, Bidentate C18 2.0, Replacement Cartridges, 2.0mm ID x 10mm Length, 2.2um, 120A. Hichrom style, in individual black cases. Cogent TYPE-C Brand. 5/PK.



Additional Info

These Guard Columns are designed to protect your 2.0 analytical columns from particulate and contaminants that will adsorb to your column. These columns are packed in Stainless Steel tubes with SS Frits. Labelled, and shipped in individual, black, Delrin cases to protect and keep them from drying out. Use with [Hichrom Style Guard Column Holders](#).

Use these guard columns with 2.1mm ID analytical columns packed with 2.2um Bidentate C18.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,091.83

Printed On 2024/05/19