


Catalog Number

45008-Y

Tubing Connector Nuts, low pressure, 1/8 inch, Polypropylene, Yellow with 1/4-28 screw threads, for use with flanged, polymeric tubing. ARE-Applied Research brand. 10/PK.



Pack of 10

Additional Info

These tubing connector nuts are used to connect 1/8th inch (3.2 mm) polymeric tubing in the low pressures zones of HPLC systems and other low to medium pressure, fluidic uses.

- Compatible with 1/4 - 28 UNF fitting ports.
- Useful up to approximately 500 psi or 34.5 bar.
- Maximum Temperature: 100°C.
- Click [HERE](#) for a convenient flanging tool.

Tubing must be flanged, and washers are required to work properly, and sold separately, click [HERE](#) for washer ordering information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$42.38

Printed On 2024/05/18