

Catalog Number

45099-002

Flanging Tool Kit with 5 Drills, Tubing Cutter, Tubing Connectors, Tubing & Plastic Carrying Case. ARE-Applied Research brand. 1 EA.



Additional Info

"Flangy the Tool" for cold flanging of Polymer or thin walled PEEK Tubing, provides an easy and smooth way to make or restore, tubing tip, flanges easy. As the flange is created by mechanical means rather than by heat, the resulting flange itself will be more long lasting and stable.

This Kit includes the Flanging Tool, set of 5 flanging drills (pins), a sample of PTFE Tubing that is 1/16" ID x 0.75" OD, six feet long, and the Clean Cut Tubing Cutting Tool, a selection of Tubing Connectors (nuts) with PP rings and the carrying case.

(5) Flanging Drills (pins)

- PEEK Tubing Drill with 0.5 mm SS Pin [45099-105](#)
- Polymer Tubing Drill with 0.8 mm Plastic Pin [45099-108](#)
- Polymer Tubing Drill with 0.8 mm Titanium Pin [45099-308](#)
- Polymer Tubing Drill with 1.3 mm Plastic Pin [45099-113](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$553.73

Printed On 2024/05/18