

Catalog Number

46000-LP

Fittings Kit for HPLC. low pressure. Includes and assortment of flange free fittings and ferrules, in an organizing box. ARE-Applied Research brand. 1 EA.



Additional Info

This convenient kit puts all the parts you need to keep the low pressure zones of an HPLC up and running. A must have for everyone that is responsible for maintenance of the instruments.

Includes:

- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Black 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Blue 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Green 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Red 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", White 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/16", Yellow 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Black 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Blue 1/4-28
- 10 Connectors, Tubing, Delrin, Flange Free, 1/8", Green 1/4-28

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$617.76

Printed On 2024/05/18