


## Catalog Number

47305-ST-10-3M

**Tubing, PEEK, 0.010 inch (0.25 mm) ID, 1/16th inch (1.6 mm) OD,**  
Ultra-T grade, blue intermittent stripes, 3 meter roll, ARE-Applied  
Research brand. 1 EA.



### Additional Info

This Ultra-T™ grade (*premium*) tubing is extruded from virgin PEEK material without any manufacturing additives. The broken blue stripes are added only to the outside of the tubing, leaving the interior all natural PEEK which eliminates leachables due to color.

- ID and OD tolerances: +/- 0.025 mm
- Wall Thickness: 0.675 mm
- Color code: blue stripes
- Pressure rating: 5600 psi (386.1 bar)

Other lengths available:

- [1 meter roll](#)
- [5 meter roll](#)
- [10 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$92.45**

Printed On 2024/05/20