


Catalog Number

47310-100-1M

Tubing, PEEK, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD,
HPLC grade, natural color, 1 meter roll, ARE-Applied Research brand. 1
EA.



Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.80 mm
- Color code: natural
- Pressure rating: 3200 psi (220.6 bar)
- *Caution Note: This tubing has very thick walls and can unwind very quickly. Handle with care to prevent injury. Safety glasses and gloves should be worn when unpacking.*

Other lengths available:

- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$40.61

Printed On 2024/05/20