


## Catalog Number

47310-100-C

**Tubing, PEEK, 1/16th inch (1.6 mm) ID, 1/8th inch (3.2 mm) OD,**  
HPLC grade, natural color, 25 meter roll, ARE-Applied Research brand. 1  
EA.



### Additional Info

This wide OD, PEEK, tubing is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 0.80 mm
- Color code: natural
- Pressure rating: 3200 psi (220.6 bar)
- Supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$914.85**

Printed On 2024/05/20