

Catalog Number

47310-75-1C

Tubing, PEEK, 0.030 inch (0.75 mm) ID, 1/8th inch (3.2 mm) OD,
HPLC grade, natural color, Your custom continuous length tubing, priced
per meter. ARE-Applied Research brand.



Additional Info

You determine the length of this wide OD, PEEK, tubing. It is extruded from virgin PEEK material without adding any additives.

- ID and OD tolerances: +/- 0.05 mm
- Wall Thickness: 1.225 mm
- Color code: natural
- Pressure rating: 5000 psi (344.7 bar)
- Lengths of 15 m or longer supplied on a convenient spool.

Other lengths available:

- [1 meter roll](#)
- [3 meter roll](#)
- [5 meter roll](#)

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/20